6 Using Decision-Making Structures

In this lesson you will learn to create different decision-making structures and be able to identify applications where using these structures can be beneficial.

Topics

- + Case Structures
- + Event-Driven Programming

Exercises

Exercise 6-1 Temperature Warnings With Error Handling

Exercise 6-2 Converting a Polling Design to an Event Structure Design


A. Case Structures

Objective: Recognize and use the basic features and functionality of Case Structures.



Activity 6-1: Case Structures Review

Figure 6-1. Case Structures Review




- Case Selector Label
- 2 Selector Terminal
- 1. What is the purpose of the Case Structure?
 - Execute one of its subdiagrams based on an input value
 - b. Repeat a section of code until a condition occurs
 - c. Execute a subdiagram a set number of times
- 2. How many of its cases does a Case Structure execute at a time?
 - a. All of them
 - b. One
- 3. What is the purpose of the Case Selector Label?
 - a. Lets you wire an input value to determine which case executes
 - b. Show the name of the current state and enable you to navigate through different cases
- 4. What is the purpose of the Selector Terminal?
 - a. Lets you wire an input value to determine which case executes
 - b. Show the name of the current state and enable you to navigate through different cases



Activity 6-1: Case Structures Review - Answers

Figure 6-2. Case Structures Review: Answers




Case Selector Label

- 2 Selector Terminal
- What is the purpose of the Case Structure?
 - a. Execute one of its subdiagrams based on an input value
 - Repeat a section of code until a condition occurs
 - Execute a subdiagram a set number of times
- 2. How many of its cases does a Case Structure execute at a time?
 - All of them a.
 - b. One
- 3. What is the purpose of the Case Selector Label?
 - Lets you wire an input value to determine which case executes
 - b. Show the name of the current state and enable you to navigate through different cases
- 4. What is the purpose of the Selector Terminal?
 - Lets you wire an input value to determine which case executes
 - b. Show the name of the current state and enable you to navigate through different cases




Demonstration: Case Structures


Right-click the Case structure to display the shortcut menu. The shortcut menu gives you options for configuring a Case structure.



Selector Terminal Data Types

You can wire a variety of data types to the Selector Terminal. The Case structure configuration changes based on the type of data that you connect.






Input and Output Tunnels

As with other types of structures, you can create multiple input and output tunnels.

- Input tunnels are available to all cases if needed.
- Output tunnels require that you define a value for each case.
- Be cautious using the Use Default If Unwired option.
 - Adds a level of complexity to the code.
 - Can complicate debugging your code.
 - The default value may not be the value that you expect.





- A tunnel that has been wired for all cases.
- A tunnel that has not been wired.
- 2 A tunnel that is marked as Use Default If Unwired. Right-click on the tunnel and select Use Default if Unwired.

Data Type Wired to Tunnel	Default Value
Numeric	0
Boolean	FALSE
String	empty ("")

Demonstration: Selector Terminal Types and Tunnels

- Create Case structures using different data type selectors.
- Create different types of output tunnels.





Exercise 6-1 Temperature Warnings With Error Handling

Goal

ů

Modify a VI to use a Case structure to make a software decision.


Scenario

You created a VI where a user inputs a temperature, a maximum temperature, and a minimum temperature. A warning string generates depending on the temperature that is less than the minimum temperature. Modify the VI to generate a different string to alert the user to the error: Upper Limit < Lower relationship of the given inputs. However, a situation could occur that causes the VI to work incorrectly. For example, the user could enter a maximum Limit. Set the Warning? indicator to TRUE to indicate the error.

Design

Modify the flowchart created for the original Temperature Warnings VI as shown in Figure 6-3.

Figure 6-3. Modified Temperature Warnings Flowchart



You must add a Case structure to the Temperature Warnings VI to execute the code if the maximum temperature is less than or equal to the minimum temperature. Otherwise, the VI does not execute the code. Instead, the VI generates a new string and the Warning? indicator is set to TRUE.




Figure 6-4. Original Temperature Warnings VI Block Diagram

Implementation

- Open Weather Warnings.lvproj in the <Exercises>\LabVIEW Core 1\Weather Warnings directory.
- Open Temperature Warnings.vi from the Project Explorer window.
- Open the block diagram and create space to add the Case structure.

To select more than one item press the <Shift> key while you select the items.

Select the Weather Data In type-defined cluster terminal, the Unbundle by Name function, and the Error In terminal.

While the objects are still selected, use the left arrow key on the keyboard to move the controls to the left.

- Press and hold the <Shift> key to move the objects in five pixel increments.




Press the <Ctrl> key and use the Positioning tool to drag out a region of the size you want to insert.

Select the Weather Data Out type-define cluster terminal, the Bundle by Name function, and the Error Out terminal.

- While the terminals are still selected, use the right arrow key on the keyboard to move the indicators to the right.
- Select the wire connecting the Weather Data In terminal and the Bundle by Name function.
- While the wire is still selected, use the up arrow key on the keyboard to move the wire upward.
- Modify the block diagram similar to that shown in Figure 6-5, Figure 6-6, and Figure 6-7. This VI is part of the temperature weather station project. 4.

Temperature Warnings VI Block Diagram-No Error, False Case Figure 6-5.



- Less?—Compares the Max Temperature and Min Temperature. Make sure the Less? function is outside the Case structure.

 Case Structure—Do not include the Weather Data In, Error In, Weather Data Out, or Error Out terminals in the Case structure because these controls - 2
- and indicators are used by both cases.
 Set True and False cases—With the True case visible, right-click the border of the Case structure and select **Make this Case False**.

 Case Structure—Wire the Error In terminal to the selector terminal to create No Error and Error cases. By default, the Case structure has True and False cases. These cases change to Error and No Error cases only after you wire Error In to the selector terminal.

Weather Data Out Error Out (7) (\mathfrak{o}) If Max Temperature setting < Min Temperature setting, output a warning ■ True ▼ No Error ▼ § <u>√</u> Error In Weather Data In

Figure 6-6. Temperature Warnings VI-No Error, True Case

- True case—If the Max Temperature is set lower than the Min Temperature, the True case executes. Click the case selector label to choose the True
 - True Constant—When the True case executes, the Temperature Warning? LED illuminates in the Weather Data Out cluster.

 String Constant—If the Max Temperature is set lower than the Min Temperature, the warning Upper Limit < Lower Limit displays on the front panel. Enter the text in the String Constant. 20
- Predict the values for Temperature Warning Text and Temperature Warning? given each set of inputs. 2

Table 6-1. Predict Values for Temperature Warnings VI

Current Temperature	Max Temperature	Min Temperature	Temperature Warning Text	Temperature Warning?
30	30	10		
25	30	10		
10	30	10		
25	20	30		

Create the Error case in the outer Case structure so this VI can be used as a subVI. 9.

Weather Data Out Figure 6-7. Temperature Warnings VI—Error Case Ė Weather Data In

7. Save the VI.

Test

- 1. Switch to the front panel of the VI.
- Test the VI by entering values from Table 6-2 in the Current Temperature, Max Temperature, and Min Temperature controls and running the VI for each set of data. 5

Table 6-2 shows the expected Temperature Warning Text and Temperature Warning? Boolean value for each set of data.

Testing Values for Temperature Warnings VI Table 6-2.

Temperature Warning?	True	False	True	True
Temperature Warning Text	Heatstroke Warning	No Warning	Freeze Warning	Upper Limit < Lower Limit
Min Temperature	10	10	10	30
Max Temperature	30	30	30	20
Current Temperature	30	25	10	25


- Do these values match the values that you predicted?
- What happens if you set the value for all three inputs to 10?
- How could you address this problem?
- Test the Error case. To use this VI as a subVI, the VI must be able to handle an error coming into the VI. Test the Error case to make sure that this VI can output the error information it receives. რ
- On the front panel, use the Operating tool to click the status Boolean indicator inside the Error In cluster so that the indicator turns red and enter 7 in the code control.
- Run the VI. The error information you entered passes through the Error case in the VI and is output in the Error Out cluster.
- Display the block diagram, select the No Error case, highlight execution, and then run the VI again to see the error pass through the Error case.
- On the front panel, right-click the border of the Error Out cluster and select Explain Error to display information about the error that was returned.
- Save and close the VI.

End of Exercise 6-1

B. Event-Driven Programming

Objective: Recognize basic features and functionality of event structures.

Demonstration: Event-Driven Scenario






Event-driven programming

Event

Method of programming where the program waits on an event to occur before executing the code written to handle that event.

An asynchronous notification that something has occurred. Events can originate from the user interface, external I/O, or other parts of the program. In this course, you will only learn about user interface events, which include mouse clicks, key presses, and value changes of a control.


Polling Versus Events






Multimedia: Event-Driven Programming

Complete the multimedia module, Event-Driven Programming, available in the <Exercises>\ LabVIEW Core 1\Multimedia\ folder to learn about programming with events.




Configuring the Event Structure

You can select which events the Event structure implements by right-clicking the border and selecting Edit Events Handled by This Case from the shortcut menu





Edit Events Dialog Box



Configured events


Event sources

Events

Notify and Filter Events

LabVIEW categorizes user interface events into two different types of events:

- Notify (Green arrow) Notify events inform you that a user action occurred. LabVIEW has already performed the default action associated with that event.
- Filter (Red arrow) Filter events allow you to validate or change the event data before LabVIEW performs the default action associated with that event. You also can discard the event entirely to prevent the change from affecting the VI.







Note A single case in the Event structure cannot handle both notify and filter events. A case can handle multiple notify events but can handle multiple filter events only if the event data items are identical for all events.

Demonstration: Configure and Use Events

Configure and use an Event structure to create a VI that responds to user interface events using event-driven programming.





Exercise 6-2 Converting a Polling Design to an Event Structure Design

Goal

ő

To convert a polling-based application to an event-based application and compare the different in performance.

Description

First you observe the behavior of a polling VI.

Next, you modify the polling VI to create a more efficient, event-driven VI and observe the changes in behavior.

Finally, you add different types of events to the VI.

Table 6-3 lists the events you will implement in the UI Event Handler VI you create.

Table 6-3. User Interface Events

Event Description	Stops the While Loop.	Displays a time stamp when you click the Time Check button.	Displays the coordinates of the front panel point you click.	Handles the event in which the user tries to close the running VI by clicking the window close button.	Produces a beep when the mouse cursor moves over the Stop button.
Event	"Stop": Value Change	"Time Check": Value Change	"Pane": Mouse Down	Panel Close?	"Stop": Mouse Enter


Observing the Polling VI Behavior

- Open and run Polling.vi.
- Open the Events.lvproj file in the <Exercises>\LabVIEW Core 1\Events directory and open the Polling VI from the project.
- 2. Examine the performance of a polling VI using the Windows Task Manager.
- ☐ Press the <Ctrl-Alt-Delete> keys and select Start Task Manager from the menu.
- ☐ Click the Performance tab in the Windows Task Manager window.
- ☐ Run the VI.

	☐ Notice how high the CPU usage is.
	☐ Stop the VI and notice how the CPU usage drops.
რ	Open the block diagram, turn on execution highlighting, and run the VI again.
4.	Notice how often the Time Check terminal sends data to the Case structure and how often the While Loop iterates.
5.	Stop the VI and turn off execution highlighting.
≥ ←	Modifying the Polling VI to Use Events Instead of Polling 1. Save Polling VI as UI Event Handler.vi so you can modify it.
	Select Open additional copy and add the copy to the project.
2	Close Polling.vi.
က်	Open the block diagram of UI Event Handler.vi and move the Stop terminal and the Time Check terminal outside the While Loop. You move these terminals into the appropriate event cases later in this exercise.
4.	Delete the Case structure and clean up any broken wires.
5.	Place an Event structure inside the While Loop between the iteration terminal and the conditional terminal.
9.	Right-click the Event structure and select Edit Events Handled by This Case from the shortcut menu.

Configure the event as shown in Figure 6-8. 7.

Figure 6-8. Configuring the "Stop": Value Change Event




Click Stop in the Event Sources panel. Click Value Change in the Events panel.

Click **OK** to close the dialog box.

Place a True constant inside the new "Stop": Value Change event and wire it to the conditional terminal of the While Loop as shown in Figure 6-9. ი ი

Figure 6-9. Event Structure with "Stop": Value Change Event




NewVal event data—Resize the event data items list so that only one item displays. Click the item and select NewVal.

Observing the Behavior of the Event-Driven VI

- 1. Run the VI.
- 2. Notice that the Iteration indicator does not increment.
- 3. Switch to the block diagram and enable execution highlighting.
- Notice that the While Loop is executing the first iteration. The Event structure is waiting for an event. 4.
- 5. Disable execution highlighting and switch back to the front panel.
- 6. Click the **Stop** button to stop the VI.
- 7. Notice that the VI stops running even though the Stop button is disconnected.
- Notice that the Stop button stays depressed even though the mechanical action is set to Latch When Released. The reason the button stays depressed is because the VI stopped running after you clicked the button. ω.
- 9. Reset the Stop button by clicking it again.

10. Drag the terminal of the Stop button into the "Stop": Value Change event as shown in Figure 6-10.

Figure 6-10. "Stop": Value Change Event with Stop Button Terminal




11. Run the VI and click the Stop button again.

12. Notice this time the VI stops and the button resets.

Programming the "Time Check": Value Change Event

- 1. Add a new event case and create a "Time Check": Value Change event as shown in Figure 6-11.
- Right-click the event structure and select Add Event Case.

Figure 6-11. Event Structure with "Time Check": Value Change Event




- In the Edit Events window, select Time Check in the Event Sources panel and Value Change in the Events panel − 0 € 4
 - Move the Time Check terminal from outside the While Loop into the "Time Check": Value Change event case. Get Date/Time In Seconds—Creates a time stamp in memory. Indicator—Displays the current time output of the Get Date/Time In Seconds function.
- Run the VI.
- Click the Time Check button to see the current time display in the current time indicator. ω.
- Display the Task Manager window and notice that CPU usage has decreased when you use events instead of polling. 4.
- Stop the VI. ъ.

Adding More Notify Events to the VI

1. Add a new event case and create a Mouse Down event as shown in Figure 6-12.

Event Structure with "Pane": Mouse Down Event Figure 6-12.




- In the Edit Events window, select Panes»Pane in the Event Sources panel and Mouse»Mouse Down in the Events panel. 7 0 0
- Coords event data—Click the event data node and select Coords»All Elements.

 Coords indicator—Right-click the output of the Coords event data item and select Create»Indicator from the shortcut menu.
- Run the VI.
- Click on different parts of the front panel. რ
- Notice that the Coords indicator displays the coordinates for each point you click.
- Notice that the other events continue to behave as before.
- Stop the VI. 4.

Adding Filter Events to the VI

1. Add a new event case and create a Panel Close? event as shown in Figure 6-13.

Event Structure with Panel Close? Event Figure 6-13.



After you add the event, in the Edit Events window, select < This VI> in the Event Sources panel and Panel Close? in the Events panel

Two Button Dialog function and Not function—Wire the T button? output to the Not function and wire the Not function to the Discard? event filter node. data node—Click the Event Data Node and select Source from the menu. Event

String constant—Wire Are you sure you want to close the window? to the message input. −284G

Yes and No string constants—Wire Yes to the T button name ("OK") input and wire No to the F button name ("Cancel") input.

- Save and run the VI.
- Click the "X" at the top-right of the window of the front panel. ω.
- Notice that clicking the No button cancels the event and returns to the VI. 4.
- Clicking the Yes button stops and closes the VI. Я.
- Stop the VI if necessary 6.

Challenge

1. If you have a sound card, add an event that produces a sound when the cursor is over the Stop button.




Tip Use Quick Drop to find the Beep VI

End of Exercise 6-2

Caveats and Recommendations


The following list describes some of the caveats and recommendations to consider when incorporating events into LabVIEW applications.

- Place only one Event structure in a loop.
- Use a Value Change event to detect value changes.
- Keep event handling code short and quick.
- Place Boolean control terminals inside an event case for latched operations to work properly.
- Avoid using an Event structure outside of a loop.
- Avoid configuring two Event structures for the same event.



Think about the VIs that you will need to develop at your job.

Will you use event-based programming to implement any of your VIs? Why or why not?



Learn More About	LabVIEW Help Topic or ni.com
Event-driven programming	Locking Front Panels
	Choosing How the Event Structure Monitors For Events
	Viewing Enqueued Events at Run Time
	Event-Driven Programming
	Events in LabVIEW
Deciding whether to use events	Caveats and Recommendations when Using Events in LabVIEW
Part of the Event structure	Event Structure
	Configuring Events Handled by the Event Structure
Determining which notifier to use	Determining Which Type of User Interface Events to Use



Activity 6-2: Lesson Review

- 1. Which of the following can NOT be used as the case selector input to a Case structure?
 - a. Error cluster
 - b. Array
 - c. Enum
 - d. String
- 2. How many events can an Event structure handle each time it executes?
 - a. As many as have occurred since the last time the event structure executed
 - b. One per configured event case
 - c. One
- 3. Which statements about event-driven programming versus polling are true?
 - a. Events execute on demand.
 - b. Event-driven programming is less CPU-intensive.
 - c. Event structures handle all events in the order the occur.
 - d. Polling may fail to detect a change.



Activity 6-2: Lesson Review - Answers

- 1. Which of the following can NOT be used as the case selector input to a Case structure?
 - a. Error cluster
 - b. Array
 - c. Enum
 - d. String
- 2. How many events can an Event structure handle each time it executes?
 - a. As many as have occurred since the last time the event structure executed
 - b. One per configured event case
 - c. One
- 3. Which statements about event-driven programming versus polling are true?
 - a. Events execute on demand.
 - b. Event-driven programming is less CPU-intensive.
 - c. Event structures handle all events in the order the occur.
 - d. Polling may fail to detect a change.